

EasyBag Classic

Isothermal bag

for the transport of insulin or adrenaline
at an optimal temperature

EasyBag Classic features

THE FREEDOM TO TRAVEL

Designed for the transport of five prefilled insulin or adrenaline pens at an optimum temperature between 16 to 25°C for 5 days

COOL WITHOUT ELECTRICITY

The EasyBag Classic becomes active in less than 40 seconds when the inner pouch is soaked in plain tap water. As long as the inner pouch is wet, your medication is kept at the right temperature

RECOMMENDED USES

Pre-filled pens of insulin or adrenaline

PATENTED TECHNOLOGY BY BASF IN GERMANY

CAN BE RE-USED HUNDREDS OF TIMES

QUICK AND EASY TO USE

IDEAL FOR ACTIVE TRAVELLERS

DISCREET AND ELEGANT

*Discover all the MedActiv solutions
for the safe transport of medication
www.medactiv.com*

Imagined and developed in France by
HealthWorks France
75008, Paris, France
contact@medactiv.com



The EasyBag Classic keeps your insulin cool for 5 days and is the perfect pouch to keep your insulin at the right temperature (between 16 and 25°C) for longer trips



Technical specifications

- Weight: 90g
- Dimensions: 210 x 160 x 35mm
- Operating Environment:
Temp: 0 °C ~ 32 °C, humidity: ≤85% RH
- Conditions for transport and storage:
Temp: -20 °C ~ 55 °C, humidity <95% RH



medactiv
transforming lives